

AIR FUMIGANT

DARMOSAN 5

- Bactericidal
- Tuberculocidal

- Compatible with Equipment Surfaces
- Formulated with Patented Hydrogen Peroxide Technology
- No Harsh Chemical Fumes

- Virucidal*
- Fungicidal

- One-Step Disinfectant Fumigant
- No Harsh Chemical Fumes
- Non-Flammable



**Dry Fog
Leaves No Residue**

**← AUTOMATIC FUMIGATING
FOGGER**

With
**Active
Oxygen**

5 kg

KEEP OUT OF REACH OF CHILDREN

READ THE LABEL BEFORE USING

STANDARD OPERATING PROCEDURE FOR AEROJET STERIZ.

- Switch off the AHU of the sterility room and inoculation room, wear the mask, gloves and goggles before starting the work.
- Loosen tank clamp by pulling finger tabs upwards.
- Lift the power head of the tank and set it down on a clean surface.
- Add fogging liquid (DARMOSAN 5) to the tank-MAXIMUM 4LITER
- Replace and re-clamp the power head of the tank
- Set the required flow control with the help of flow control valve
- Lock the valve with rotating red nut for constant flow of liquid solution out of the fogger.
- Place it on Roto fog turn table
- Plug the fogger machine and Roto fog turntable in mains
- Place it at appropriate place in the room.
- Switch on the mains so that the machine generates the fog.
- After completion of fogging treatment (45mins for 1000cu.ft), remove the tank and add sterile distilled water into the tank and operate for 2-3mins to ensure cleaning of internal parts. Switch off the machine and switch off the mains.
- Clean the tank.

APPLICATIONS AREAS

- Pharmaceutical manufacturing areas
- Laboratories
- Operation theatres, ICU, clinics etc
- Disinfection measures for diseases and epidemic control
- Hatcheries, poultry shades, animal houses.
- Medical devices
- Food processing plants and dairies.
- Hotel rooms.
- Green houses.
- Closed units etc.

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. **FIRST AID: If in eyes:** Hold eye open and rinse slowly and gently with water for 15–20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. Call the doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor or going for treatment.

STORAGE AND DISPOSAL: Store this product in a cool, dry area, away from direct sunlight and heat to avoid deterioration. Nonrefillable container. Do not reuse or refill empty container. Recycle empty container or discard in trash.

ACTIVE INGREDIENT:

Hydrogen Peroxide,
Colloidal silver (Stabilizer)
Dilutant

Batch No.:

Best before:

MICROBES SANITIZATION

- Acinetobacter baumannii
- Campylobacter jejuni
- Carbapenem Resistant Klebsiella pneumoniae
- Enterobacter aerogenes
- Escherichia coli (E. coli)
- Escherichia coli New Delhi Metallo-Beta Lactamase-1 (NDM-1)
- Escherichia coli O157:H7
- Extended Spectrum Beta Lactamase producing Escherichia coli (ESBL producing E. coli)
- Extended Spectrum Beta Lactamase producing Klebsiella pneumoniae (ESBL producing Klebsiella pneumoniae)
- Klebsiella oxytoca
- Klebsiella pneumoniae New Delhi Metallo-BetaLactamase-1 (NDM-1)
- Multidrug Resistant Enterococcus faecium
- Multidrug Resistant Klebsiella pneumoniae
- Pseudomonas aeruginosa
- Salmonella enterica
- Stenotrophomonas maltophilia
- Penicillin-Resistant Streptococcus pneumoniae
- Streptococcus pyogenes
- Vancomycin Resistant Enterococcus faecalis (VRE)

BACTERIA: 1 minute contact time

- Staphylococcus epidermidis (Coagulase-negative staphylococci) (CoNS)
- Community Acquired Methicillin Resistant Staphylococcus aureus (CA-MRSA)
- Linezolid Resistant Staphylococcus aureus
- Methicillin Resistant Staphylococcus aureus (MRSA) (ATCC 33592)
- Staphylococcus aureus
- Staphylococcus aureus (MRSA) (ATCC 33591)
- Methicillin Resistant Staphylococcus aureus (MRSA)

BACTERIA: 4 minute contact time

- Mycobacterium bovis (TB)

*VIRUSES: 30 second contact time

- *Avian Influenza A virus (H5N1)
- *Cytomegalovirus
- *Hepatitis B virus
- *Hepatitis C virus
- *Herpes simplex virus type 1
- *Herpes simplex virus type 2
- *Human Coronavirus
- *Human Immunodeficiency Virus type 1 (HIV-1)
- *Influenza A virus (H3N2)
- *Influenza A virus (2009 H1N1)
- *Influenza B virus
- *Respiratory syncytial virus (RSV)

*VIRUSES: 1 minute contact time

- *Norovirus
- *Rhinovirus (type 37)
- *Rotavirus

MOLD, FUNGI: 3 minute contact time

- Candida albicans
- Trichophyton mentagrophytes

MANUFACTURED BY:



OSHO CHEMICAL INDUSTRIES LTD.

P. O. Box 49916 00100 Nbi, Kenya
Pilot lines: [+254] 0711 045000 | 0732 167000 | 0726 880 480 | 0733 880 480
SMS: 20560 | Fax: [+254] 0711 045600 E-mail: oshochem@oshochem.com
Web: www.oshochem.com